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[Menu] → [Setting] → [Assay Parameters] → [PT THS] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol							
PT THS	Rinse at the Time of Primary											
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse											
Fibrinogen	Sample Dispensing /Dilution											
Thrombine Time	Dilution Setting											
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8		Edit	Add					
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List	Edit						
F VIII	Sample	50	0	-	-	Delete						
F IX												
BCPIg												
	Factor-Deficient Plasma					Add						
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List	Edit		Delete					
							Add					
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit					
	PT THS	100	190	Weak	-	Clean I	Delete					
	Measurement Condition											
	Measurement Time (Main)	180 sec	Measurement Time (Sub)	180 sec	Seq.	Reagent	Lot management	Stable				
					1	OVB	<input type="checkbox"/>	0				
	Measurement Time (Main)	Measurement Time (Sub)			2	PT THS	<input checked="" type="checkbox"/>	0				
	<input type="checkbox"/>	Mixing during measurement										
	STEP 1	0 sec –	0	0 rpm								
	STEP 2	0 sec –	0	0 rpm								
	STEP 3	0 sec –	180 sec	0 rpm								
	Export							Add	Load	Print	Save	Close

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[Menu] → [Setting] → [Assay Parameters] → [APTT PSL] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol			
PT THS	Rinse at the Time of Primary							
APTT PSL	Aspiration <input type="checkbox"/>		Pre-Rinse <input type="checkbox"/>	Post-Rinse				
Fibrinogen	Sample Dispensing							
Thrombine Time	/Dilution							
AT III	Dilution Setting						Edit	
BCPC	Buffer:	None	Dilution Ratio	1/1				
F VIII	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List	Edit		
F IX	Sample	50	0	-	-	Delete		
BCPIg								
	Factor-Deficient Plasma							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List	Edit			
					Delete			
	Add							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	
	APTT PSL	50	60	Weak	-	-	Delete	
	CaCl2	50	190	Weak	-	-		
	Measurement Condition							
	Measurement Time (Main)	180 sec	Measurement Time (Sub)	180 sec	Seq.	Reagent	Lot management	Stable
					1	APTT PSL	<input checked="" type="checkbox"/>	0
					2	CaCl2	<input type="checkbox"/>	0
	Measurement Time (Main)	Measurement Time (Sub)						
	<input type="checkbox"/>	Mixing during measurement						
	STEP 1	0 sec –	0	0 rpm				
	STEP 2	0 sec –	0	0 rpm				
	STEP 3	0 sec –	180 sec	0 rpm				
	Export							
	Add							
	Load							
	Print							
	Save							
	Close							

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[Menu] → [Setting] → [Assay Parameters] → [Fbg] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol			
PT THS	Rinse at the Time of Primary							
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse							
Fibrinogen	Sample Dispensing /Dilution							
Thrombine Time	Dilution Setting							
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 2/1, 3/1		Edit	Add	
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List	Edit		
FVIII	Sample	10	90	-	-	Delete		
FIX								
BCPIg								
	Factor-Deficient Plasma							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List	Edit			
	Delete							
	Add							
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	
	Fbg	50	190	Weak	-	Clean I	Delete	
	Measurement Condition							
	Measurement Time (Main)	100 sec	Measurement Time (Sub)	100 sec	Seq.	Reagent	Lot management	Stable
	Measurement Time (Main)	Measurement Time (Sub)		1	OVB	<input type="checkbox"/>	0	
	<input type="checkbox"/> Mixing during measurement			2	Fbg	<input checked="" type="checkbox"/>	0	
	STEP 1	0 sec -	0	0 rpm				
	STEP 2	0 sec -	0	0 rpm				
	STEP 3	0 sec -	100 sec	0 rpm				
	Export							
	Add							
	Load							
	Print							
	Save							
	Close							

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[Menu] → [Setting] → [Assay Parameters] → [ТТ] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol			
PT THS	Rinse at the Time of Primary							
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse							
Fibrinogen	Sample Dispensing /Dilution							
Thrombine Time	Dilution Setting							
AT III	Buffer:	OVB	Dilution Ratio	1/1	Edit	Add		
BCPC	Aspirated Sample	Sample Aspiration Vol [µl]	Buffer Volume [µl]	Pre-rinse List	Post-rinse List	Edit		
FVIII	Sample	50	0	-	-	Delete		
FIX								
BCPIg								
Factor-Deficient Plasma								Add
Aspirated Reagent		Reagent Aspiration Vol [µl]	Pre-rinse List	Post-rinse List				Edit
								Delete
								Add
Aspirated Reagent		Reagent Aspiration Vol [µl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	
Test Thr		100	70	Weak	-	Clean I	Delete	
Measurement Condition								
Measurement Time (Main)		160 sec	Measurement Time (Sub)	160 sec	Seq.	Reagent	Lot management	Stable
Measurement Time (Main)		Measurement Time (Sub)						
<input type="checkbox"/>		Mixing during measurement						
STEP 1	0 sec –	0	0 rpm					
STEP 2	0 sec –	0	0 rpm					
STEP 3	0 sec –	160 sec	0 rpm					
					1	OVB	<input type="checkbox"/>	0
					2	Test Thr	<input checked="" type="checkbox"/>	0
Export		Add	Load	Print	Save	Close		

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[Menu] → [Setting] → [Assay Parameters] → [AT3] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol			
PT THS	Rinse at the Time of Primary							
APTT PSL	Aspiration	<input type="checkbox"/>	Pre-Rinse	<input type="checkbox"/>	Post-Rinse			
Fibrinogen	Sample Dispensing							
Thrombine Time	/Dilution							
	Dilution Setting							
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 3/2, 0/1		Edit	Add	
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List		Edit	
FVIII	Sample	21	147	-	-		Delete	
FIX	Dilution	24	0					
BCPIg								
	Factor-Deficient Plasma					Add		
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List			Edit	
							Delete	
							Add	
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	
	AT3Thro	175	40	Weak	Clean I	Clean I	Delete	
	AT3Subs	33	230	Strong				
	Measurement Condition							
	Measurement Time (Main)	60 sec	Measurement Time (Sub)	60 sec	Seq.	Reagent	Lot management	Stable
					1	OVB	<input type="checkbox"/>	0
	Measurement Time (Main)		Measurement Time (Sub)		2	AT3Thro	<input checked="" type="checkbox"/>	0
	<input type="checkbox"/>	Mixing during measurement			3	AT3Subs	<input checked="" type="checkbox"/>	0
	STEP 1	0 sec –	0	0 rpm				
	STEP 2	0 sec –	0	0 rpm				
	STEP 3	0 sec –	60 sec	0 rpm				
	Export	Add	Load	Print	Save	Close		

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[Menu] → [Setting] → [Assay Parameters] → [BCPIg] → [Test Protocol]

Assay Group	Basic	Re-analysis	Reflex	QC	Test Protocol				
PT THS	Rinse at the Time of Primary								
APTT PSL	Aspiration <input type="checkbox"/> Pre-Rinse <input type="checkbox"/> Post-Rinse								
Fibrinogen	Sample Dispensing /Dilution								
Thrombine Time	Dilution Setting								
AT III	Buffer:	OVB	Dilution Ratio	1/1, 1/2, 1/4, 1/8, 3/2, 0/1			Edit	Add	
BCPC	Aspirated Sample	Sample Aspiration Vol [μl]	Buffer Volume [μl]	Pre-rinse List	Post-rinse List			Edit	
FVIII	Sample	21	147	-	-			Delete	
FIX	Dilution	28	0						
BCPIg									
	Factor-Deficient Plasma								
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pre-rinse List	Post-rinse List			Edit	Add	
	PIReag						Delete		
	Add								
	Aspirated Reagent	Reagent Aspiration Vol [μl]	Pipette at [sec]	Mixing Level	Pre-rinse List	Post-Rinse List	Edit	Delete	
	PIReag	175	40	Middle	Clean I	Clean I			
	PISubs	35	350	Strong					
	Measurement Condition								
	Measurement Time (Main)	120 sec	Measurement Time (Sub)	120 sec	Seq.	Reagent	Lot management	Stable	
					1	OVB	<input type="checkbox"/>	0	
	Measurement Time (Main)	Measurement Time (Sub)			2	PIReag	<input checked="" type="checkbox"/>	0	
	<input type="checkbox"/>	Mixing during measurement				3	PISubs	<input checked="" type="checkbox"/>	0
	STEP 1	0 sec –	0	0 rpm					
	STEP 2	0 sec –	0	0 rpm					
	STEP 3	0 sec –	120 sec	0 rpm					
	Export	Add	Load	Print	Save	Close			

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