

## Номера реагентов для штрикодов на IL ACL TOP

<b>Код реагента</b>	<b>Наименование реагента</b>	<b>Тип флакона</b>	<b>Код теста</b>
701	Renamplastin	2	401
702	APTT Renam	2	402, 410
703	APTT CaCL <sub>2</sub> Renam	2	402, 410
704	FIB C Renam	2	403, 404, 405
705	Thrombin Renam	2	406
706	AT Thrombin Renam	2	407
707	AT Substrate Renam	2	407
708	PLG Streptokinase Re	2	408
709	PLG Substrate Ren	2	408
710	Prot C Activator Ren	2	409
711	Prot C Substrate Ren	2	409
712	Prot C Buffer Ren	2	409
713	FVIIIDefRen	2	410
714	Cal Plasma Ren	2	401,403,404, 405,407,408, 409,410
715	Plasma N	2	401,403,404, 405,407,408, 409,410
716	Plasma P	2	401,403,404, 405,407,408, 409,410